



FLOW CYTOMETRY BUFFERS FOR CYTOKINE RESEARCH

Lymphocyte Stimulation Cocktail

Fixation / Permeabilization

Cell Staining Buffer

Enhance the quality of your fluorescent immunostaining experiments with these pre-optimized flow cytometry buffers from SouthernBiotech.

These buffers are particularly suited to cytokine research. The Lymphocyte Stimulation Cocktail ensures optimum accumulation of intracellular cytokines by increasing protein synthesis and cellular retention. Following cellular stimulation - the Fixation/Permeabilization Buffer Kit prepares cells for intracellular staining in 40 minutes. Finally, the Cell Staining Buffer may be used to stain both intact and permeabilized cells.

Buffer Type	Cat. No.		Size
Lymphocyte Stimulation Cocktail	0210-01		200µl
Fixation / Permeabilization Buffer	0215-01	Fixation / Permeabilization Buffer	50mL
		Permeabilization Wash Buffer	50mL
Cell Staining Buffer	0225-01S		100mL
	0225-01M		250mL
	0225-01L		500mL