


## MOUNTING SOLUTIONS AND SLIDES FOR IHC

**Fluoromount-G™**

**Dapi-Fluoromount-G™**

**Gelatin Subbed Slides**



Improve your immunohistochemistry work flow with these mounting solutions from SouthernBiotech.

Process your samples on our ready-to-use slides which are suitable for IHC embedding and detection techniques. The slides come in two convenient sizes for low and high throughput use.

Staining with antibodies for fluorescence microscopy will benefit from the two mounting media we offer:

**Fluoromount-G™** and **Dapi-Fluoromount-G™** are both aqueous media that offer convenient slide mounting, with a semi-permanent seal for long-term storage of slide preparations.

**Dapi-Fluoromount-G™** additionally carries out nuclear staining with DAPI (455nm)

Cytokine Research Reagents	Cat. No.		Size
Fluoromount-G™	0100-01	Clear liquid mounting medium	25 mL
Dapi-Fluoromount-G™	0100-20	Clear liquid mounting medium	20 mL
Gelatin Subbed Slides	SLD01-BX	Premium slides subbed with Gelatin	72 slides/box
	SLD01-CS	Premium slides subbed with Gelatin	10 boxes/case