



Welcome to Nordic Labs' first company newsletter. Our aim is to keep you informed of our top selling products, new products, and news that we think may be of interest to you.

## PRODUCT FOCUS

### ANTI-MONKEY ANTIBODIES

Antibodies are available specific for an IgG, IgA, and IgM common determinant as well as the individual isotypes. Antibodies can be chosen for their Fc, Fab, or heavy and light chain binding. In addition to whole serum and purified antibody, the following formats are available: Biotin, FITC, HRP and TRITC.

All our anti-monkey antibodies are free of preservatives and supplied lyophilized.

#### Adsorbed Secondary antibody:

The majority of Nordic Lab's antibodies have been immunoaffinity adsorbed using insolubilized antigens to eliminate:

1. Antibodies cross reacting with other serum proteins (especially removal of Ig/Fab in Fc specific antibodies)
2. Prevent the presence of excess adsorbent protein or immune complexes in the antiserum.

#### Fab Fragment of affinity purified secondary antibodies:

Generated by papain digestion of whole IgG. Ideal for direct staining of fixed cells, IHC, and ELISA.

#### Fc Fragment of affinity purified secondary antibodies:

Ideal for direct Immunofluorescence staining of cytoplasmic Ig, electrophoresis, and double radial immunodiffusion.

### Anti-Monkey Immunoglobulins

|                  | Host   | Serum | Unconjugated | Biotin | FITC | HRP | TRITC |
|------------------|--------|-------|--------------|--------|------|-----|-------|
| IgA(Fc specific) | Goat   | ✓     | ✓            | ✓      | ✓    | ✓   | ✓     |
|                  | Rabbit | ✓     |              |        | ✓    |     |       |
| IgG (H&L)        | Goat   | ✓     | ✓            | ✓      | ✓    | ✓   |       |
|                  | Rabbit | ✓     | ✓            | ✓      | ✓    | ✓   |       |
| IgM(Fc Specific) | Goat   | ✓     | ✓            | ✓      | ✓    | ✓   | ✓     |
|                  | Rabbit | ✓     |              |        |      |     |       |
| IgG Fab          | Rabbit | ✓     | ✓            | ✓      | ✓    | ✓   | ✓     |
| IgG(Fc)          | Goat   | ✓     | ✓            | ✓      | ✓    | ✓   | ✓     |
|                  | Rabbit | ✓     | ✓            | ✓      | ✓    | ✓   |       |

### Other Monkey Specificities

|   | Host   | Serum | Unconjugated | Biotin | FITC | HRP | TRITC |
|---|--------|-------|--------------|--------|------|-----|-------|
| Anti-albumin (IgG)                      | Rabbit | ✓     |              |        |      |     |       |
| Anti-DNP-albumin (polyclonal)           | Rabbit | ✓     |              |        |      |     |       |
| Anti-C3c (IgG)                          | Goat   | ✓     | ✓            | ✓      | ✓    | ✓   | ✓     |
| Anti-lactoferrin (IgG)                  | Goat   | ✓     | ✓            | ✓      | ✓    | ✓   | ✓     |
| Anti-Secretory component (free & bound) | Goat   | ✓     | ✓            | ✓      | ✓    | ✓   | ✓     |
| Serum Proteins                          | Rabbit | ✓     |              |        |      |     |       |

## TOP SELLING PRODUCTS THIS QUARTER

| Product Code            | Description                    |
|-------------------------|--------------------------------|
| 1. MAHu/IgG4/MAB (1 mg) | Mouse Anti Human IgG4          |
| 2. GAM/IgG(H+L)/FITC    | Goat Anti Mouse IgG(H&L):FITC  |
| 3. MAHu/IgG4/PO         | Mouse Anti Human IgG4:HRPO     |
| 4. RAG/IgG(H+L)/FITC    | Rabbit Anti Goat IgG(H&L):FITC |
| 5. MAHu/IgG1/Bio        | Mouse Anti Human IgG1:Biotin   |
| 6. GARa/IgG(H+L)/FITC   | Goat Anti Rat IgG(H&L):FITC    |
| 7. MAHu/IgGc/Bio        | Mouse Anti Human IgG:Biotin    |
| 8. RAM/Ig/FITC          | Rabbit Anti Mouse Ig:FITC      |
| 9. GAHu/IgG(Fc)/PO      | Goat Anti Human IgG(Fc):HRPO   |
| 10. SwAM/IgG(H+L)/FITC  | Swine Anti Mouse IgG(H&L):FITC |

## THE NORDIC LABS QUARTERLY PROMOTION

**BUY ANY MONKEY ANTIBODY  
AND GET A  
SECOND MONKEY ANTIBODY  
AT 50 % DISCOUNT**

*\* If you wish to participate in this offer, please contact Tim Bernard, who will issue you with a promotional code:  
**tim@pivotalscientific.com***

*(Offer valid until the end of November 2011)*